

Question

The determination of glucose content in a protein sample is usually performed using spectrophotometry. Weigh 12.2 g of the sample, dissolve it in water, add saturated sodium chloride solution, filter, and dilute to 250 mL to obtain the test solution X . Prepare glucose standard solutions of different concentrations, take 5.00 mL of each solution, add 4.00 mL of Fehling's reagent, dilute to 10 mL with water, heat in a boiling water bath for 30 minutes, centrifuge, take the supernatant, and measure its absorbance at 660 nm. Treat the test solution X in the same way. The experimental data obtained are shown in the following table:

Concentration $c(\text{g} \cdot \text{L}^{-1})$	Absorbance A
0.150	0.660
0.250	0.568
0.350	0.466
0.450	0.383
0.550	0.284
0.650	0.188
X	0.348

In addition to spectrophotometry, we can also use titration to determine the glucose content. Pipette 10.00 mL of Fehling's reagent, add 10 mL of water and 5 mL of potassium ferrocyanide solution, add two glass beads, heat quickly to boiling, and use methylene blue as an indicator. Titrate rapidly with a $0.500 \text{ g} \cdot \text{L}^{-1}$ glucose standard solution until the solution becomes colorless, consuming a volume of 15.00 mL. Titrate with the test solution X as the titrant using the same method, consuming a volume of 15.80 mL.

Please combine the aforementioned analysis process, analyze the role of each step, and select the matching options.

- A.** All other options are incorrect.
- B.** The statement "Saturated sodium chloride is added during sample processing to salt out proteins" is incorrect.
- C.** The statement "In spectrophotometry, absorbance should have a linear relationship with concentration within this range" is incorrect.
- D.** The statement "In titration, the purpose of adding glass beads is to prevent bumping" is incorrect.
- E.** The statement "In titration, the role of potassium ferrocyanide is as a reducing agent" is incorrect.
- F.** The statement "When measuring a sample using titration, if the endpoint time is long, the measurement results will generally be lower" is correct.
- G.** The statement "During sample processing, the purpose of removing proteins is because they can also be oxidized by Fehling's reagent" is correct.
- H.** The statement "the described spectrophotometric method can be used for the analysis of glycoproteins" is correct.
- I.** The statement "The glucose content in the original sample measured by titration is 0.475%" is correct.
- J.** The statement "The glucose concentration in sample X measured by spectrophotometry is $0.988 \text{ g} \cdot \text{L}^{-1}$ " is correct.

Answer

Correct Answer: **F**

Detailed Explanation

First, salting out is used to remove proteins, but glycoproteins will also be salted out, and the sugar in them is not reducing sugar, so it cannot be simply treated in this way.

CHECKPOINT

0.5 PTS

Salting out is to remove proteins

CHECKPOINT

0.5 PTS

This method cannot analyze glycoproteins

Proteins will exhibit a color reaction with the biuret reagent of the Fehling's reagent, so they need to be separated, but this is not a redox reaction.

CHECKPOINT

1 PTS

Proteins exhibit a color reaction with the biuret reagent of the Fehling's reagent

Next, spectrophotometry is analyzed using the standard curve method, which requires being within the linear range.

CHECKPOINT

1 PTS

Standard curve method analysis requires a linear relationship between absorbance and concentration

Next, through the analysis of known samples, the slope of the linear standard curve between concentration (g/L) and absorbance is -0.941 , and the intercept is -0.803 .

CHECKPOINT

0.5 PTS

Standard curve slope is -0.941 , intercept is -0.803

Substituting the absorbance data of X , the concentration of sample X measured by spectrophotometry is $0.482 \text{ g} \cdot \text{L}^{-1}$, and the glucose content in the original sample is 0.988% .

CHECKPOINT

1 PTS

The concentration of sample X measured by spectrophotometry is $0.482 \text{ g} \cdot \text{L}^{-1}$, and the glucose content in the original sample is 0.988%

In the titration method, since there is a boiling operation, adding glass beads is a common method to prevent bumping.

CHECKPOINT

0.5 PTS

Adding glass beads is to prevent bumping

Potassium ferrocyanide in this system is used to mask the color of copper ions for observing the endpoint.

CHECKPOINT

1 PTS

Potassium ferrocyanide is used to mask color

The titration process needs to be carried out quickly, because air will also oxidize the reduced product of methylene blue back, so the added glucose in the titration is too much, and the result will be smaller.

CHECKPOINT

1 PTS

Titration should be performed quickly to avoid air oxidation of the catalyst, resulting in a smaller titration result

Then, from the titration results, the concentration is inversely proportional to the consumption volume, and the concentration of sample X is $0.475 \text{ g} \cdot \text{L}^{-1}$, and the glucose content in the original sample is 0.973%.

CHECKPOINT

1 PTS

The concentration of sample X measured by titration is $0.475 \text{ g} \cdot \text{L}^{-1}$, and the glucose content in the original sample is 0.973%

